

B Safe

High efficiency sanitizing spray system




COMPRESSORS
medical

B/Safe is MGF's new high efficiency sanitizing spray system.

B/Safe utilizes a high performance turbine for a superior efficient dry-spraying of disinfection solutions.

Manufactured in Italy and with decades of presence on the market, **B/Safe** is your **simple and reliable solution for ambient sanitizing!**

TANK

NOZZLE FOR OPTIMAL NEBULIZATION



ADJUSTMENTS FOR LIQUID DISPENSING

QUALITY & PERFORMANCE

The B/Safe turbine provides a high airflow spray-mist at low pressure with the following benefits:

- **Constant air flow** for optimal spraying, without pressure fluctuation
- **High spray range:** up to 5 feet
- **High efficiency:** only 17 oz. of sanitizing liquid for the treatment of 3,2 cubic feet!
- **Dry spraying:** due to the calibrated nozzle orifice and warm turbine air, there is no need to dry the surfaces after the treatment.

The professional air gun is made with high quality materials:

- 0,5 mm nozzle for the optimal nebulization
- Simple settings for the optimal use of the disinfection liquid
- An adjustable cap allowing for **different operating modes according to the types of surfaces to be treated**
- The 17 oz. plastic gravity fed gun permits the efficient use of disinfection solutions.

TECHNICAL FEATURES

Codice prodotto:

- **TB-055-S:** turbine
- **TB-055-S-L:** version equipped with a 1 lt bottle of **Atom Pur** sanitizing liquid

Power: **550W (115V-60Hz)**

Quantity of air delivered: **64 cfm**

Pressure: **16,4 psi**

Dimension: **11,8x8,7x12,2 in**

Weight: **10,5 lb**

Operating temperature: **-10°C/+40°C**

Noise level: **75 dB**

EASY TO USE

B/Safe is handy, light and ready for use!

- **Less than 11 lb** for best handling
- **The electric cable and the air pipe are both 9,9 feet**, for the easiest operating
- Comes complete with a **shoulder belt** and an **accessories holster**, where to carry the airbrush when not in use
- **Adjustment free turbine:** switch it on and it is **ready for use!** Simply adjust the spray flow level on the airgun

DISINFECTION LIQUID SELECTION

Never use **B/Safe** with alcohol-based chemicals or other dangerous or highly flammable volatile liquid components.

MGF recommends the use of the sanitizing liquid **Atom Pur**, reference codes:

- **TB-055-S-L-1:** 1 lt bottle
- **TB-055-S-L-10:** 10 lt bottle

Otherwise, MGF recommends the use of hydrogen peroxide, as recommended by the World Health Organization and most of the sanitary disinfection protocols.

For further details, please contact your dealer.



MGF Compressors s.r.l.
Via Pascoli 15
20081 Abbiategrasso (MI) ITALY
Tel +39 02 90 19 180
Fax +39 02 90 19 273
Email export@mgfcompressors.it
Web www.mgfcompressors.com